


Usage and Effect

The ear canal naturally massages ear wax & other debris into the outer ear using tiny hairs called "Cilia" that move when you talk or chew. Occasionally wax can back up in the auditory canal, this can commonly be caused by using ear plugs, hearing aids or the use of cotton buds.

flip-tips® are designed to prevent blockages by working with the ears natural, self-cleaning process by gathering ear wax & skin particles on its soft, cleaning blades - leaving a sheilding layer of fat for acid protection against pathogenic germs.

Deeper intrusions into the auditory canal are reserved only for specialist, medical care personnel.



Recommended Usage of flip-tips®:

1. Insert flip-tips® slowly into your ear canal, ensure you hold it from the rubberized handle
2. As flip-tips® move into the ear being to slow rotate clockwise
3. Once you have fully rotated flip-tip clockwise, begin to rotate anti-clockwise
4. Slowly remove flip-tips® and clean by rising under the tap

flip-TIPS!

- Slow rotating improves the effect
- Wet the blade for removing very dry cerumen
- Repeat every 2-3 days as described
- We recommend replacing your flip-tips® every 3-6 months